



GLOBAL
MEDICAL ENDOSCOPY



GME-LS10
Electrosurgical Energy Platform
Vessel Sealing Energy

Product Features

INTELLIGENT TISSUEFFECT™ PLUS TECHNOLOGY

Applies appropriate energy to targeted tissue, for an effective seal that will withstand a minimal 3x the normal systolic blood pressure.

FAST SEALING

Seals in 2-4 seconds.

SIMPLE CONTROL PANEL

Single power on/off button, single LigaSure™ instrument plug-in port.

SMART CONNECTOR

Recognizes which LigaSure™ instrument is being used and automatically configures energy settings for quick, consistent results.

COMPACT & LIGHTWEIGHT

Uses minimal space and is easily portable.

Size: 30.0 x 10.5 x 37.7 cm.

Weight: 5.5 kg.



Green - Tseal



Flashing - Sealing in process



Purple - ST Tseal



Yellow background - Information

Intelligent bloodless cutting



- Isolate tissue
- Scrape tissue
- Coagulated tissue
- Vessels Sealing
- Widely used in a variety of surgeries
- Compatible With Covidien Energy Platform

Product Code	Head Type	Seal Plate Length	Shaft Length	Cut Length	Shaft Diameter
CY1923A		20.55mm	23cm	18.65mm	5mm
CY1937A		20.55mm	37cm	18.65mm	5mm
CY1944A		20.55mm	44cm	18.65mm	5mm
CY1923B		20.55mm	23cm	18.65mm	5mm
CY1937B		20.55mm	37cm	18.65mm	5mm
CY1944B		20.55mm	44cm	18.65mm	5mm



- Blunt tip vessel sealer and divider
- Low temperature profile Coagulated tissue
- Target tissue with thermal spread < 2mm
- Vessels/Rectum Sealing
- Compatible With Covidien Energy Platform

Product Code	Head Type	Seal Plate Length	Shaft Length	Cut Length	Shaft Diameter
CY1020		20.00mm	20cm	20.00mm	10mm
CY1037		20.00mm	37cm	20.00mm	10mm



- Blunt tip vessel sealer and divider
- Grasp, seal, and cut tissue independently
- Minimize the need for multiple activations
- Vessels/Rectum Sealing
- Compatible With Covidien Energy Platform

Product Code	Head Type	Seal Plate Length	Shaft Length	Cut Length	Shaft Diameter
CY4418		36.00mm	18cm	34.00mm	The Nonstick Coating



- Blunt tip vessel sealer and divider
- Scrape tissue
- Coagulated tissue
- Vessels/Tumor/Tissue Sealing
- Compatible With Covidien Energy Platform

Product Code	Head Type	Seal Plate Length	Shaft Length	Cut Length	Shaft Diameter
CY1212		17.30mm	19cm	15.20mm	For use in open procedures Such as Ear, Nose and Throat (ENT)

Intelligent bloodless cutting



GLOBAL

MEDICAL ENDOSCOPY

Endoscopy Equipment Specialists



**Entrega a
disponibilidad**



**Certificado de
garantía**



**Equipos médicos
de calidad**

Contactanos

 909-563-8510

23891 Via Fabricante
Suite 609
Mission Viejo,
CA 92691 USA



 www.globalmedicalendoscopy.com