

LT-300 HD

High Definition Digital Video Colposcope

Features

- SONY digital camera
- 2.38 Mega Pixels
- Mini LCD viewing screen
- Adjustable LED lights
- Progressive green filter
- Auto and manual focus
- Handheld remote
- Foot pedal for hands-free image capture



Applications




Certifications

FDA **CE** ISO 13485:2012 **Made in USA**



Mounting Options



Vertical Stand



Swing Arm

Specifications

Camera	SD: CMOS	HD: SONY Exmor CMOS
Mega Pixels	SD: 2MP	HD: 2.38MP
Magnification	SD: Optical: 1-20X, Digital: 1-30X HD: Optical: 1-30X, Digital: 1-120X	
Visual Angle	360° camera rotation	
Working Distance	SD: 7.1-14.9in	HD: 5.1-15.7in
Video Output	TV: HDMI	Computer: USB 3.0
Green Filter	Progressive green filter offers a wider range of contrast	
Focus	Automatic and manual focus	
Image Freeze	Image freeze on screen available with remote, foot pedal, or on camera	
Remote	Remote minimizes motion while adjusting the image	
Battery	Can work up to 4.5 hours using battery power	
Foot Pedal	Foot pedal available for hands free picture and video	



Laptop



Monitors



Desktop



EHR & EMR Server



GLOBAL
MEDICAL ENDOSCOPY

