



IPC™ System

INTEGRATED POWER CONSOLE



Innovation *that Transforms*

Inspired by Surgeons and Nurses in the Field

The Integrated Power Console is the **ONLY** powered ENT system to offer:

- The widest range of application-specific products
- IntelliFlow™ remote control for irrigation
- 0.5 mm/minute irrigation flow rates for working in the facial recess
- Integrated facial nerve monitoring with the Visao® Otologic Drill and Stim Bur Guard*
- High-speed, application-specific curved burs for otologic procedures
- A complete anterior skull base surgery solution

*Nerve monitoring capability provided by our NIM-Response® 2.0 or NIM-Neuro® 2.0 nerve monitors.

Integrated Power Console (IPC™) System *Innovation*

The new Integrated Power Console delivers innovation that transforms your O.R. Its design was inspired by feedback from you, the surgeons and nurses in the field. Our IPC™ system is the only powered ENT surgery console with the most expansive range of application-specific products. From inferior turbinoplasties to complex skull base cases, the IPC platform is your ENT surgery system.

*Integrated facial nerve monitoring with the Visao® Otologic Drill and Stim Bur Guard**

Easy-to-use touchscreen menu for each handpiece

IntelliFlow™ remote control for irrigation

Maximum visualization for anterior skull base

Connect multiple handpieces at once



*Nerve monitoring capability provided by our NIM-Response® 2.0 or NIM-Neuro® 2.0 nerve monitors.

...n that Transforms

Complete Integration

Power all of your favorite Medtronic ENT handpieces

- Straightshot® M4 Microdebrider
- Endo-Scrub® 2 Lens Cleaning Sheaths
- Midas Rex® Legend EHS Stylus™ Drill
- Skeeter® Otologic Drill
- Visao® High-Speed Otologic Drill



Flexible Irrigation

On-drill irrigation reinvented

- Disposable IntelliFlow™ remote allows precise flow control from sterile field
- 0.5 mm/minute irrigation flow rates for working in the facial recess
- Snap-on cassette for easier set-up
- Inner diameter of tubing enables continuous flow, even at low flow rates



Multifunction Footpedal

Ergonomic control at your feet

- Three buttons provide maximum control
- Built-in footrests and ergonomic angle for comfort
- Hoop locks firmly in place, then collapses easily for storage
- Non-slip design and compact size is ideal for your O.R.



Inferior Turbinate Blade

Powered inferior turbinoplasty

- Offers significantly better long-term results compared to electro-surgical techniques^{1,2}
- Very rarely requires repeat procedures, unlike radiofrequency^{1,3}
- Reduces the risk of complications, such as pain, crusting and osteitis²



Anterior Skull Base Surgery

The complete solution for superior visualization

- Visao® High-Speed Otologic Drill
- Curved burs and blades allow better visibility and access
- Integrated suction keeps surgical site clear

The IPC platform is the ONLY powered ENT system that provides the widest range of application-specific products.

Stim Bur Guard for the Visao® Otologic Drill

Integrated facial nerve monitoring

- Provides 1-3 mm advance warning of facial nerve proximity
- Helps reduce the risk of patient injury for added peace-of-mind
- Works with our NIM® intraoperative nerve monitor*



*For use with Medtronic's NIM-Response® 2.0 and NIM-Neuro® 2.0

Ordering Information

IPCSYSKIT

IPC™ System

includes:

Product	Description	Qty
1898001	Integrated Power Console	1
1898430	Integrated Power Console Multifunction Footpedal	1
1897821	Power Cord, 6 Meter, IEC 320, 115V	1
1898851	Integrated Power Console Manual	1
1897510	Basket	1

IPCES2SYSKIT

IPC™ System and Endo-Scrub® 2 Accessories

includes:

Product	Description	Qty
1898001	Integrated Power Console	1
1898430	Integrated Power Console Multifunction Footpedal	1
1897821	Power Cord, 6 Meter, IEC 320, 115V	1
1898851	Integrated Power Console Manual	1
1897510	Basket	1
1852000	Footswitch Endo-Scrub 2	1
1991015	Finger Switch Endo-Scrub 2	1

M4KIT1

Straightshot® M4 Microdebrider and IPC™ Accessories

includes:

Product	Description	Qty
1898200T	XPS® Straightshot® M4 High-Speed IGS Microdebrider	1
1898400	XPS Straightshot M4 Sterilization Tray	1
