


OLYMPUS

Your Vision, Our Future

HYSTEROSCOPES OES 4000


Slim outer Diameter, excellent Image Quality.



- Complete line of OES 4000 stainless steel sheaths
- Perfect flow characteristics
- Scratch-proof sapphire glass front window
- Ultra-wide field of view

3 mm **standard** system

4 mm **advanced** system



Perfect flow characteristics







Complete line of OES 4000 stainless steel sheaths











The 3 mm **standard** system provides you with the hysteroscope of choice, an optical quality close to the 4 mm telescope, great ease of use. For extensive therapeutic procedures in the operation room we recommend our 4 mm **advanced** system. The use of high value materials and a unique design secure longterm durability and Olympus superior autoclavability.








Standard System

3 mm Telescope	
A4674A	Telescope, 3 mm, 30° direction of view, wide angle, autoclavable  
 <p>Delivery includes an instrument tray and an autoclavable protection sheath.</p>	
WA03300A	Light guide cable, 2.8 mm, plug type 

Advanced System

4 mm Telescope	
A22005A	Telescope, 4 mm, 30° direction of view, wide angle, autoclavable  
 <p>Delivery includes an instrument tray.</p>	
WA03300A	Light guide cable, 2.8 mm, plug type 

Sheaths	
A4770	Sheath, single flow, 4 mm 
A4771	Sheath*, continuous flow, 4.5 mm, 3 Fr. channel 
A4772	Sheath*, continuous flow, 5.5 mm, 5 Fr. channel 
A4773	Sheath*, continuous flow, 6.5 mm, 7 Fr. channel 

Sheaths and Working Inserts	
A4760	Sheath, single flow, 5 mm 
A4761	Sheath*, continuous flow, 5.5 mm, 3 Fr. channel 
A47610A	Sheath*, continuous flow, 6.5 mm, 5 Fr. channel 
A4767	Obturator, standard, for A4760, A4761, A47610A 
A4763	Sheath, 8 mm 
A4766	Obturator, 8 mm, optical, for A4763 
A4764	Working insert*, with ramp, 7 Fr. channel 

* For sealing caps and surgical instruments (HYS+), see our catalog or contact your local representative.

Specifications, design, and accessories are subject to change without any notice or obligation on the part of the manufacturer.