




GLOBAL

MEDICAL ENDOSCOPY



Endoscopy Equipment Specialists



FICHA TÉCNICA

Mini Nephroscope

Mini Nephroscope




Mini Nephroscope

Model	Diameter	Working Instrument	Working Length	Match Product
JY2256A	12Fr	6Fr	225mm	8968.421
JY2506A	12Fr	6Fr	250mm	8968.405
JY2706A	12Fr	6Fr	270mm	Custom-made
JY2756A	8.5/12Fr	6Fr	275mm	Custom-made

Nephroscope Sheath

Model	Diameter	Working Length	Match Product
JY2256B	15Fr	205mm	JY2256A
JY2506B	15Fr	230mm	JY2506A
JY2706B	15Fr	250mm	JY2706A
JY2756B	15Fr	255mm	JY2756A





A handle

R5410G3 Rigid grasping forceps
5Fr/410mm



B handle



R6410T Rigid saw-tooth forceps
6Fr/410mm



R6410T Rigid foreign body forceps
6Fr/410mm



R6410T Rigid grasping forceps
6Fr/410mm



R6410T Rigid biopsy forceps
6Fr/410mm





C0620P Puncture cannula
6Fr/206mm



C1113D Dilating cannula
11Fr/135mm



T1035S Suction tube
10.5Fr/350mm



D0822 Dilator
8-22.5Fr



R9580D Rigid dilator rod
9Fr/580mm



F9580D Flexible dilator rod
9Fr/580mm



R1340 Cleaning rod
1.5mm/340mm





GLOBAL

MEDICAL ENDOSCOPY

Endoscopy Equipment Specialists



**Entrega a
disponibilidad**



**Certificado de
garantía**



**Equipos médicos
de calidad**

Contactanos

 909-563-8510

23891 Via Fabricante

Suite 609

Mission Viejo,

CA 92691 USA



 www.globalmedicalendoscopy.com